

Image analysis to quantify Dvl2_mEos condensates

Christian Scheeder¹ and Antionia Schubert¹

¹Division Signalling and Functional Genomics, German Cancer Research Center (DKFZ) Heidelberg

December 10, 2021

Contents

1	Load dependencies	2
2	Preamble	2
3	Load the exmaple images	2
4	Perform image analysis	5
5	Session info	11

1 Load dependencies

```
knitr::opts_chunk$set(echo = TRUE)
knitr::opts_chunk$set(fig.align="center")
knitr::opts_chunk$set(dev="png")
library(EImage)
library(tidyverse)
## -- Attaching packages ----- tidyverse 1.3.1 --
## v ggplot2 3.3.5      v purrr 0.3.4
## v tibble 3.1.6       v dplyr 1.0.7
## v tidyr 1.1.4        v stringr 1.4.0
## v readr 2.1.1        v forcats 0.5.1
## -- Conflicts ----- tidyverse_conflicts() --
## x dplyr::combine()   masks EImage::combine()
## x dplyr::filter()    masks stats::filter()
## x dplyr::lag()       masks stats::lag()
## x purrr::transpose() masks EImage::transpose()
library(here)
## here() starts at /Users/c.scheeder/Desktop/remotes/Supp_Schubert_2021/condensate_quantification
```

2 Preamble

Script to document the image analysis carried out to quantify condensates. In brief, The image analysis was carried out as follows: In a first step cells were identified by segmentation. Nuclei were identified based on a Hoechst staining (DAPI channel) and cell bodies based on a DyLight Phalloidin staining (Cy5 channel). After segmentation condensates inside cells were identified based on intensity-based thresholding (FITC channel) and counted. Condensate counts and cell counts (number of segmented cells) were saved for each image. Four different cell lines were compared: - HEK cells with a DVL1-mEOS tag - HEK cells with a DVL1-mEOS tag and a EVI k.o. (clone 1) - HEK cells with a DVL1-mEOS tag and a EVI k.o. (clone 3) - HEK cells with a DVL1-mEOS tag and a FZD1 k.o. (clone 1) - HEK wild-type The image analysis shown here is exemplary carried out for two cell lines (DVL1-mEOS and DVL1-mEOS_EVIk.o._clone1).

3 Load the exmaple images

```
source_dir_DVL1_mEOS <-
  here("./raw_data/example_images/HC10984t_HC1098_E01_Evi1_1_2020.10.21.23.32.00")


source_dir_DVL1_Evi1 <-
  here("./raw_data/example_images/HC10984t_HC1098_E01_mEos_2_2020.10.21.15.56.44")

images_DVL1_mEOS <- EImage::combine(
  readImage(
    file.path(source_dir_DVL1_mEOS,
              list.files(
                source_dir_DVL1_mEOS)[grepl("DAPI",
```

Image analysis to quantify Dvl2_mEos condensates

```
list.files(source_dir_DVL1_mEOS))]]),  
  readImage(  
    file.path(source_dir_DVL1_mEOS,  
              list.files(  
                source_dir_DVL1_mEOS)[grepl("Cy5",  
                                              list.files(source_dir_DVL1_mEOS))]]),  
    readImage(  
      file.path(source_dir_DVL1_mEOS,  
                list.files(  
                  source_dir_DVL1_mEOS)[grepl("FITC",  
                                                list.files(source_dir_DVL1_mEOS))]])  
    )  
  ## Warning in readTIFF(x, all = all, ...): TIFFReadDirectory: Unknown field with  
  ## tag 65400 (0xff78) encountered  
  
  ## Warning in readTIFF(x, all = all, ...): TIFFReadDirectory: Unknown field with  
  ## tag 65400 (0xff78) encountered  
  
  ## Warning in readTIFF(x, all = all, ...): TIFFReadDirectory: Unknown field with  
  ## tag 65400 (0xff78) encountered  
  
  display(normalize(images_DVL1_mEOS),  
           method = "raster",  
           all=T)
```

Image analysis to quantify Dvl2_mEos condensates




```
images_DVL1_Evi1 <- EBImage::combine(  
  readImage(  
    file.path(source_dir_DVL1_Evi1,  
      list.files(  
        source_dir_DVL1_Evi1)[grepl("DAPI",  
                                     list.files(source_dir_DVL1_Evi1))]),  
  readImage(  
    file.path(source_dir_DVL1_Evi1,  
      list.files(  
        source_dir_DVL1_Evi1)[grepl("Cy5",  
                                     list.files(source_dir_DVL1_Evi1))]),  
  readImage(  
    file.path(source_dir_DVL1_Evi1,  
      list.files(  
        source_dir_DVL1_Evi1)[grepl("FITC",  
                                     list.files(source_dir_DVL1_Evi1))])  
  )  
## Warning in readTIFF(x, all = all, ...): TIFFReadDirectory: Unknown field with  
## tag 65400 (0xff78) encountered  
  
## Warning in readTIFF(x, all = all, ...): TIFFReadDirectory: Unknown field with
```

Image analysis to quantify Dvl2_mEos codensates

```
## tag 65400 (0xff78) encountered

## Warning in readTIFF(x, all = all, ...): TIFFReadDirectory: Unknown field with
## tag 65400 (0xff78) encountered

display(normalize(images_DVL1_Evil),
        method = "raster",
        all=T)
```



4 Perform image analysis

```
# blur the images
images_DVL1_mEOS_smooth <- images_DVL1_mEOS

images_DVL1_mEOS_smooth[, , 1] <- gblur(images_DVL1_mEOS_smooth[, , 1],
                                       radius = 51,
                                       sigma = 1)
```

Image analysis to quantify Dvl2_mEos condensates

```
images_DVL1_mEOS_smooth[, , 2] <- gblur(images_DVL1_mEOS_smooth[, , 2],
                                         radius = 51,
                                         sigma = 4)

images_DVL1_mEOS_smooth[, , 3] <- gblur(images_DVL1_mEOS_smooth[, , 3],
                                         radius = 51,
                                         sigma = 4)

images_DVL1_Evil_smooth <- images_DVL1_Evil

images_DVL1_Evil_smooth[, , 1] <- gblur(images_DVL1_Evil[, , 1],
                                         radius = 51,
                                         sigma = 1)

images_DVL1_Evil_smooth[, , 2] <- gblur(images_DVL1_Evil[, , 2],
                                         radius = 51,
                                         sigma = 4)

images_DVL1_Evil_smooth[, , 3] <- gblur(images_DVL1_Evil[, , 3],
                                         radius = 51,
                                         sigma = 4)

# segment the nuclei
segmentNuclei <- function(inputImage){
  nucleusTresh = thresh(inputImage,
                        w = 20, h = 20,
                        offset = 0.004)

  nucleusTresh = fillHull(opening(nucleusTresh,
                                  kern=makeBrush(9, shape="disc")))

  nucleusFill = fillHull(thresh(inputImage,
                                w = 30, h = 30,
                                offset = 0.001))


  nucleusRegions = propagate(inputImage,
                             seed=bwlabel(nucleusTresh),
                             mask=nucleusFill)

  return(nucleusRegions)
}

nuclei_DVL1_mEOS <- segmentNuclei(inputImage=images_DVL1_mEOS_smooth[, , 1])
nuclei_DVL1_Evil <- segmentNuclei(inputImage=images_DVL1_Evil_smooth[, , 1])
```

Image analysis to quantify Dvl2_mEos codensates


```
display(  
  EBImage::combine(  
    colorLabels(nuclei_DVL1_mEos),  
    colorLabels(nuclei_DVL1_Evi1)  
  ),  
  all=T,  
  method = "raster"  
)
```



```
# Get outlines of cell bodies  
getCellBodies <- function(inputImage,nucleiMask){  
  cytoplasmThresh = thresh(inputImage,  
    w = 28,  
    h = 28,  
    offset = 0.001)  
  
  cytoplasmOpening = opening(cytoplasmThresh,  
    kern=makeBrush(9,shape="disc"))  
  
  cytoplasmOpening2 = opening(inputImage > 0.025)  
  
  nucleusRegions2 <- nucleiMask  
  nucleusRegions2[nucleiMask != 0] <- 1  
  
  cytoplasmCombined = cytoplasmOpening | cytoplasmOpening2 | nucleusRegions2  
  
  storage.mode(cytoplasmCombined) = "integer"  
  
  cytoplasmRegions = propagate(x = inputImage,  
    seeds = nucleiMask,
```

Image analysis to quantify Dvl2_mEos codensates

```
        lambda=1e-04,  
        mask=cytoplasmCombined)  
  
    cytoplasmRegions = fillHull(cytoplasmRegions)  
  
    return(cytoplasmRegions)  
}  
  
cellBodies_DVL1_mEOS <- getCellBodies(inputImage=images_DVL1_mEOS_smooth[, , 2],  
                                     nucleiMask = nuclei_DVL1_mEOS)  
  
cellBodies_DVL1_Evi1 <- getCellBodies(inputImage=images_DVL1_Evi1_smooth[, , 2],  
                                       nucleiMask = nuclei_DVL1_Evi1)  
  
display(  
  EBImage::combine(  
    colorLabels(cellBodies_DVL1_mEOS),  
    colorLabels(cellBodies_DVL1_mEOS)  
  ),  
  all=T,  
  method = "raster"  
)
```



```
# display the fluorescence image as normalized and pseudo-colored RGB  
# with segmentation borders painted  
ImgColor_DVL1_mEOS = rgbImage(2*normalize(images_DVL1_mEOS_smooth[, , 2]),  
                             2*normalize(images_DVL1_mEOS_smooth[, , 3]),  
                             2*normalize(images_DVL1_mEOS_smooth[, , 1]))  
  
ImgOut_DVL1_mEOS = paintObjects(cellBodies_DVL1_mEOS,  
                                paintObjects(nuclei_DVL1_mEOS,  
                                              ImgColor_DVL1_mEOS,  
                                              col='yellow'),
```



Image analysis to quantify Dvl2_mEos condensates

```
col='white')

ImgColor_DVL1_Evi1 = rgbImage(2*normalize(images_DVL1_Evi1_smooth[, , 2]),
                              2*normalize(images_DVL1_Evi1_smooth[, , 3]),
                              2*normalize(images_DVL1_Evi1_smooth[, , 1]))

ImgOut_DVL_Evi1 = paintObjects(cellBodies_DVL1_Evi1,
                              paintObjects(nuclei_DVL1_Evi1,
                                             ImgColor_DVL1_Evi1,
                                             col='yellow'),
                              col='white')

display(
  EBImage::combine(
    ImgOut_DVL1_mEOS,
    ImgOut_DVL_Evi1
  ),
  all=T,
  method = "raster"
)
```



```
# Detect mEos condensates in images by adaptive thresholding
getCondensates <- function(inputImage, cellBodyMask){
  condensates_adaptTresh <- thresh(inputImage,
                                   w = 20, h = 20,
                                   offset = 0.006)

  condensates_adaptTresh[cellBodyMask == 0] <- 0

  condensates_adaptTresh_out <- colorLabels(bwlabel(condensates_adaptTresh))


  return(condensates_adaptTresh_out)
}
```

Image analysis to quantify Dvl2_mEos condensates

```
DVL1_mEOS_conTresh <- getCondensates(inputImage = images_DVL1_mEOS_smooth[, , 3],
                                     cellBodyMask = cellBodies_DVL1_mEOS)

DVL1_Evi1_conTresh <- getCondensates(inputImage = images_DVL1_Evi1_smooth[, , 3],
                                     cellBodyMask = cellBodies_DVL1_Evi1)

display(
  EBImage::combine(
    DVL1_mEOS_conTresh,
    DVL1_Evi1_conTresh
  ),
  method = 'raster',
  all=T
)
```



```
# the number of condensates corresponds to the number of objects in
# the binary image with labels

print(paste0(
  "Number of detected condensates in the image with DVL1-mEOS cells: ",
  length(unique(bwlabel(DVL1_mEOS_conTresh)[bwlabel(DVL1_mEOS_conTresh)!=0])),
  " in ",
  length(unique(cellBodies_DVL1_mEOS[cellBodies_DVL1_mEOS!=0])),
  " cells."
))
## [1] "Number of detected condensates in the image with DVL1-mEOS cells: 20 in 973 cells."

print(paste0(
  "Number of detected condensates in the image with DVL1-mEOS Evi1 k.o. cells: ",
  length(unique(bwlabel(DVL1_Evi1_conTresh)[bwlabel(DVL1_Evi1_conTresh)!=0])),
  " in ",
  length(unique(cellBodies_DVL1_Evi1[cellBodies_DVL1_Evi1!=0])),
  " cells."
))
## [1] "Number of detected condensates in the image with DVL1-mEOS Evi1 k.o. cells: 15 in 973 cells."
```

Image analysis to quantify Dvl2_mEos condensates

```
length(unique(cellBodies_DVL1_Evi1[cellBodies_DVL1_Evi1!=0])),  
  " cells."  
)  
)  
## [1] "Number of detected condensates in the image with DVL1-mEOS Evi1 k.o. cells: 11 in 779 cells."
```

5 Session info

```
sessionInfo()  
## R version 4.1.2 (2021-11-01)  
## Platform: x86_64-apple-darwin17.0 (64-bit)  
## Running under: macOS Big Sur 10.16  
##  
## Matrix products: default  
## BLAS: /Library/Frameworks/R.framework/Versions/4.1/Resources/lib/libRblas.0.dylib  
## LAPACK: /Library/Frameworks/R.framework/Versions/4.1/Resources/lib/libRlapack.dylib  
##  
## locale:  
## [1] en_US.UTF-8/en_US.UTF-8/en_US.UTF-8/C/en_US.UTF-8/en_US.UTF-8  
##  
## attached base packages:  
## [1] stats graphics grDevices utils datasets methods base  
##  
## other attached packages:  
## [1] here_1.0.1 forcats_0.5.1 stringr_1.4.0 dplyr_1.0.7  
## [5] purrr_0.3.4 readr_2.1.1 tidyr_1.1.4 tibble_3.1.6  
## [9] ggplot2_3.3.5 tidyverse_1.3.1 EBImage_4.36.0 BiocStyle_2.22.0  
##  
## loaded via a namespace (and not attached):  
## [1] Rcpp_1.0.7 locfit_1.5-9.4 lubridate_1.8.0  
## [4] lattice_0.20-45 fftwtools_0.9-11 png_0.1-7  
## [7] rprojroot_2.0.2 assertthat_0.2.1 digest_0.6.29  
## [10] utf8_1.2.2 R6_2.5.1 tiff_0.1-10  
## [13] cellranger_1.1.0 backports_1.4.0 reprex_2.0.1  
## [16] evaluate_0.14 httr_1.4.2 pillar_1.6.4  
## [19] rlang_0.4.12 readxl_1.3.1 rstudioapi_0.13  
## [22] magick_2.7.3 rmarkdown_2.11 htmlwidgets_1.5.4  
## [25] RCurl_1.98-1.5 munsell_0.5.0 broom_0.7.10  
## [28] compiler_4.1.2 modelr_0.1.8 xfun_0.28  
## [31] pkgconfig_2.0.3 BiocGenerics_0.40.0 htmltools_0.5.2  
## [34] tidyselect_1.1.1 bookdown_0.24 fansi_0.5.0  
## [37] withr_2.4.3 crayon_1.4.2 tzdb_0.2.0  
## [40] dbplyr_2.1.1 bitops_1.0-7 grid_4.1.2  
## [43] jsonlite_1.7.2 gtable_0.3.0 lifecycle_1.0.1  
## [46] DBI_1.1.1 magrittr_2.0.1 scales_1.1.1  
## [49] cli_3.1.0 stringi_1.7.6 fs_1.5.1  
## [52] xml2_1.3.3 ellipsis_0.3.2 generics_0.1.1  
## [55] vctrs_0.3.8 tools_4.1.2 glue_1.5.1  
## [58] hms_1.1.1 jpeg_0.1-9 abind_1.4-5
```

Image analysis to quantify Dvl2_mEos condensates

```
## [61] fastmap_1.1.0      yaml_2.2.1           colorspace_2.0-2  
## [64] BiocManager_1.30.16 rvest_1.0.2          knitr_1.36  
## [67] haven_2.4.3
```